

## FLOW SORT X-RAY DIAMOND RECOVERY MACHINE TRACER TEST RECORD

SINGLE  
STAGE

DATE: YY MM DD				TIME: HH MM				FIRST STAGE	
CUSTOMER:				SORTER MODEL:				SECOND STAGE	
SORTER S/N:				CUSTOMER REF:				NO FILTER	
HR COUNT				TESTED CONDUCTED BY:				K42 FILTER	
KV		MA		RESULTS CHECKED BY:				K45 FILTER	
SELECTOR SWITCH		FEED WATER PRESSURE	kPa	MAIN WATER PRESSURE	kPa	X-RAY WATER PRESSURE	kPa		
RIGHT	SENSITIVITY SETTING	µA READING		EJECTION COUNT 1st STAGE LEFT:		EJECTION COUNT 1st STAGE RIGHT:		NUMBER OF TRACERS IN 1st STAGE CONCENTRATE BIN:	
	NUMBER OF TRACERS INPUT 1st STAGE RIGHT:		TRACERS REPORTED TO TAILINGS:						
LEFT	SENSITIVITY SETTING	µA READING		EJECTION COUNT 1st STAGE LEFT:		EJECTION COUNT 1st STAGE RIGHT:		1st STAGE CONCENTRATE BIN NUMBER	
	NUMBER OF TRACERS INPUT 1st STAGE LEFT:		TRACERS REPORTED TO TAILINGS:						

COMMENTS

## FLOW SORT X-RAY DIAMOND RECOVERY MACHINE TRACER TEST RECORD

SINGLE  
STAGE

DATE: YY MM DD				TIME: HH MM				FIRST STAGE	
CUSTOMER:				SORTER MODEL:				SECOND STAGE	
SORTER S/N:				CUSTOMER REF:				NO FILTER	
HR COUNT				TESTED CONDUCTED BY:				K42 FILTER	
KV		MA		RESULTS CHECKED BY:				K45 FILTER	
SELECTOR SWITCH		FEED WATER PRESSURE	kPa	MAIN WATER PRESSURE	kPa	X-RAY WATER PRESSURE	kPa		
RIGHT	SENSITIVITY SETTING	µA READING		EJECTION COUNT 1st STAGE LEFT:		EJECTION COUNT 1st STAGE RIGHT:		NUMBER OF TRACERS IN 1st STAGE CONCENTRATE BIN:	
	NUMBER OF TRACERS INPUT 1st STAGE RIGHT:		TRACERS REPORTED TO TAILINGS:						
LEFT	SENSITIVITY SETTING	µA READING		EJECTION COUNT 1st STAGE LEFT:		EJECTION COUNT 1st STAGE RIGHT:		1st STAGE CONCENTRATE BIN NUMBER	
	NUMBER OF TRACERS INPUT 1st STAGE LEFT:		TRACERS REPORTED TO TAILINGS:						

COMMENTS